



Product Catalogue  
 Medication Safety

Comprehensive solutions  
 for preparing, storing  
 and distributing medicines

## Improving patient safety

The WIEGAND® system offers comprehensive and flexible solutions to increase medication and, ultimately, patient safety. Parallel to providing an improvement in the quality of care, the simple, modular system components are applicable to inpatient and ambulatory healthcare environments. Our systems reduce the work load of nurses and pharmacy assistants, resulting in a reduction in medication errors while substantially decreasing costs.



### Sustainable know-how

WIEGAND® has been studying medication processes for over forty years. Through close collaboration with care specialists, pharmacists and physicians across Europe, we continually develop our concepts to find new and practical solutions for our customers' ever changing requirements.

WIEGAND®'s concepts increase efficiency and improve the working process of nurses, thus leading to a decrease in medication errors.

The WIEGAND® products are durable and have an outstanding cost-benefit ratio, which effectively facilitates the lowering of health care costs.

[www.puntogeriatrico.com](http://www.puntogeriatrico.com)





## Reducing liability risks

The frequency of medication errors occurring between writing a prescription to the administration of medicines is alarmingly high. For years, various statistics from the USA and Europe refer to this harmful and costly problem. Only by integrating these factors into a health care institution's risk management process, the reduction of medication errors will become an achievable strategic goal.

### Medication errors are costly


According to different estimates, 60-80% of all errors in the health care sector can be attributed to human failure. The consequences of medication errors are disturbing. In an extensive report by D.W. Bates and colleagues, the following extrapolations were made:

- Estimated annual costs of Adverse Drug Events (ADE) for a 700-bed hospital: \$ 5.6 million
- ADEs were recorded in 6.5% of all patients
- Overall costs: \$ 2,595 per ADE
- ADEs lead to 2.2 additional days of hospitalisation
- Preventable ADEs (28% of all ADE) annual cost: \$ 2.8 million
- Costs for each preventable ADE: \$ 4,685
- Preventable ADEs result in additional 4.6 days of hospitalisation
- 57% of all ADEs were judged significant (30% serious, 12% life-threatening and 1% fatal)

«These estimates do not include costs of injuries to patients, malpractice costs, or the costs of less serious medication errors or admissions related to ADEs. These results suggest that hospitals can justify devoting additional resources to develop systems that reduce the number of preventable ADEs not only to improve patient care but also to reduce ADE-related expenses.»


Source: Cost of adverse drug events; D.W. Bates et. al., JAMA, Jan 22/29, 1997 – Vol277, No. 4

MEDMARX® (USA) analyzed 21,084 cases with ADEs and highlighted the consequences of the harmful cases.



### Allocation of medication errors

The following data from MEDMARX® (USA) are based on 218,347 reported cases and the frequency of the different sources of error was analyzed.



Source: Hicks, R.W., Santell, J.P., Cousins, D.D., and Williams, R.L. (2004). MEDMARXSM 5th Anniversary Data Report: A Chartbook of 2003 Findings and Trends 1999–2003. Rockville, MD: USP Center for the Advancement of Patient Safety.

The Swiss Foundation for Patient Safety carried out a survey in 2006 on patient safety in Swiss hospitals. The six most pertinent sources of medication errors (number of entries x weighting) are listed below:

1. Error while preparing medicines
2. Administration of medicines to the wrong patient
3. Wrong dosage (calculation error in potencies of 10)
4. Bad legibility of drug prescriptions
5. Faulty/missing information of drug therapy with/after discharge
6. Transfer errors when copying prescriptions

Source: Stiftung für Patientensicherheit: Problemfelder (Hot-Spots) in der Patientensicherheit; Schweizerische Ärztezeitung 2008;89: 24

## Day medication prepared according to administration times

In order to guarantee a smoothly running day on the ward, the medication must be well planned and prepared in advance. Based on correct care documentation, all drugs for a whole day are put into the appropriate compartments of the medicine dispensers according to the prescribed administration times. Control of the prepared medication and their administration times is thus secured.



### Avoiding errors

- Check the patient's identity
- Ensure that dosage calculations are checked independently by another health care professional before the drug is administered
- Make sure that prescription, drug, and patient match
- Make certain that the medication is given at the correct time
- Minimise interruptions during drug rounds

Source: Medication errors, DJP Williams, J R Coll Physicians Edinb 2007; 37:343–346



# Reducing errors through the four-eyes principle

Drug therapies (medications) for all patients should be ideally prepared during a calm period. This task requires the highest level of concentration and should never take place at a patient's bedside. Before administering medicines to a patient, drugs prepared in a medication dispenser should be controlled by a second person to aid in reducing medication errors.

Medikament (MediDispenser)	07:00	08:00	12:00	18:00	21:00	Mo	Di	Mi	Do	Fr	Sa	So	Informationen	Bestand
ACTOS 300 Tabs Dosisangabe			1										Tablet, rund, weiss, nicht teilbar Ch: 110432 (1.6.05)	100.0 => 09.08.10
EUTHYROX 125 mg Tabs Dosisangabe	1												Tablet, rund, weiss, teilbar, nüchtern einnehmen Ch: 1468325 (1.6.04)	82.0 => 20.09.10
ISMN 40 retard Heumann Kap. Dosisangabe			1										Tabletten, ablang, rot-braun, nicht teilbar Ch: 240514 (1.1.08)	63.0 => 02.08.10
SILMARTIT Kaps. Dosisangabe		1	1										Wirkstoffe: Oliban, braun, nicht teilbar Ch: 9910082 (1.4.05)	28.0 => 05.07.10
METOPROLOL 200 retard Dosisangabe				1									Retard-Tablet, rund, weiss, verstellbar Ch: 000432 (1.6.07)	60.0 => 04.08.10
MAGNESIUMACID Sachet Dosisangabe					1								Pulver, in Wasser auflösen Ch: 000320 (1.1.08)	38.0 => 15.06.10
SUPRADYN Tabs. Dosisangabe			1										Brausetabletten, teilbar, in Wasser auflösen Ch: 38200002 (1.8.09)	13.0 => 19.06.10
METHYFORMIN AL 500 Dosisangabe				1									Farmakoballen, oval, weiss, halbbalbar Ch: 980102 (1.8.07)	45.0 => 21.07.10

**Wichtig**  
 Euthyrox eine halbe Stunde vor Frühstück einnehmen

**Verabreichte von:**  
 Name: \_\_\_\_\_ Unterschrift: \_\_\_\_\_

**Zusätzliche Medikamente:**  
 PREMERAN Gabe (Dosisangabe) y x

**Bedarfsmedikamente:**  
 Paracetamol 1000 mg (OP & 20 Tabs) max. Dosis: 4000mg

**Informationen:**  
 (Tabletten, nur ausserhalb anwenden) Ch: 200102 (1.2.06) 0  
 (Tabletten, oval, weiss, halbbalbar) Ch: 9150322 (1.6.08) 0

**Bestand:**  
 gerichtet: 07.06.10  
 gerichtet von: B. Smeets  
 gutig: 07.06.10 - 14.06.10  
 Lowenapotheke  
 Bahnstrasse 8  
 8000 Zurich

**Ruth Donati**  
 Geburtsdatum: 12.03.1926  
 Haus Sonnehof  
 Seestrasse 23  
 8002 Zurich  
 Station 3B  
 Zimmer 204  
 Arzt  
 Dr. Pierre Lamaze



## WIEGAND® MediDispenser

- The sliding lid opens in one direction only
- This secures the administration sequence
- 4 to 12 flexibly adjustable compartments (1 for each administration)
- Hygienic and economical (semi disposable)
- Tight sliding lid



## Adaptable medication processes

Medication processes differ substantially in terms of the logistical framework in which they are applied. A nursing home, operating with centrally or even externally prepared medications has different workflows to the medical or surgical department of a hospital. The more acute the treatments are, the shorter are the processes and operating times and thus the higher are the corresponding needs for quality in care and the costs of an incorrect treatment.



### Decentralised ward- and satellite pharmacies – cost-efficient, organised and space-saving

- Provide an excellent overview and fast access
- Low stocks and good availability due to Kanban principle (always 2 units/OP)
- Application of “first in - first out” principle of storage, provides an effective control of expiry dates and lots of the drugs
- The central pharmacy steers the quantity of supplies and determines the permitted pharmaceutical products (cost control)
- Individual drug costs are charged at time of consumption
- Fast turnover of drugs without temporary storage in BIN System



WIEGAND® Medication Prepare Station for 40 patients  
www.puntogeriatrico.com



# Improving efficiency

When using the Day System for the medication process, the nurses prepare the prescribed drugs – after the doctor’s visits – for the next 24 hours. Preparing for a 7-day period (Week System) significantly optimises operational flows in retirement- and nursing homes. To keep patients in a «normal condition», a long-term therapy with a regular, stabilizing drug therapy is applied. Therefore medications are preferably prepared for a whole week.

## Professional WIEGAND® Medication Prepare Stations

Empirical research illustrates that when using a professionally furnished Medication Prepare Station, medicines for 20 patients can be prepared in approximately one hour. Decentralised ward- and satellite pharmacies with a Medication Prepare Station are set up to manage between 40 and 100 patients.

WIEGAND® Distribution Tray  
 8 patients - 24 hours



WIEGAND® Week Tray  
 1 patient - 7 days



## BIN System combined with WIEGAND® Medication Distribution System

When using BIN Systems, a patient’s drugs are kept in individual containers. Although these systems need additional storage space, they are quite commonly used and can be easily combined with the WIEGAND® Medication Distribution System in order to be able to prepare medications in advance. In this way, patient safety can be substantially increased.



### WIEGAND® medication concepts

- Enhance existing systems, in order to increase patient safety
- Cost efficient investment to reduce the expense of medication errors
- Replace inefficient and outdated systems to reduce liability risks
- Sustainable, modular system leading to greater flexibility in application and investment security

## MediDispenser | the «small» ones

**Article No.** 01/09


**Description** MediDispenser 01

**Details**

- MediDispenser No. 01 with Sliding-Lid No. 09
- also available as imprinted Week-Set
- 165x50x20 mm

**Accessories**

510.14-8xx Label with dispensing periods  
 510.18 Label, extra-small  
 291 Week-Tray 01  
 93T-01 WeekBox 01  
 93-01 Week-Pouch 01  
 295 Sliding-Lid-Opener



**Article No.** 11/19

**Description** MediDispenser 11

**Details**

- MediDispenser No. 11 with Sliding-Lid No. 19
- 225x50x20 mm

**Accessories**

510.14-8xx Label with dispensing periods  
 510.13 Label for patients name and medication  
 510.13Wxx Label for medication with week days  
 510.13A4 Label "No. 13" for imprinting with MS-Word-form  
 591 Week-Tray 11  
 93T-11 WeekBox 11  
 91 Week-Container 11



**Article No.** 131/139


**Description** MediDispenser 131

**Details**

- MediDispenser No. 131 with Sliding-Lid No. 139
- 3 Divisions No. 132, divisible in max. 9 compartments
- 188x42x37 mm

**Accessories**

510.14-8xx Label with dispensing periods  
 510.18 Label, extra-small  
 510.17xx Label for medication  
 17SD Label holder for Label 510.17xx  
 1191 Tray for MediDispenser 131



## MediDispenser | the «medium» ones

**Article No.** 41/49


**Description** MediDispenser 41

**Details**

- MediDispenser No. 41 with Sliding-Lid No. 49-69-349
- also available as imprinted Week-Set
- 228x81x20 mm

**Accessories**

510.44xx Label with week days  
 510.46 Label, short strip  
 991 Week-Tray 41-61-341  
 193TQ-41-61-341.1bl WeekBox 41-61-341, blue design  
 191 Week-Container 41-61-341





**Article No.** 61/69

**Description** MediDispenser 61

**Details** - MediDispenser No. 61  
with Sliding-Lid No. 49-69-349  
- also available as imprinted Week-Set  
- 228x81x20 mm

**Accessories**

510.44xx Label with week days  
510.46 Label, short strip  
510.64-7 Label for medication  
66.x Inlay-Dishes, coloured  
66.Spez Singe-Use-Inlay-Dish  
67 Inlay-Blister with 4 compartments  
991 Week-Tray 41-61-341  
193TQ-41-61-341.1bl WeekBox 41-61-341, blue design  
191 Week-Container 41-61-341



**Article No.** 66.x

**Description** Inlay-Dishes, coloured

**Details** - for MediDispenser No. 61/69  
- x to be ordered in different colours:  
No. 66.1 = Inlay-Dish, transparent-white  
No. 66.2 = Inlay-Dish, transparent-yellow  
No. 66.3 = Inlay-Dish, transparent-blue  
No. 66.4 = Inlay-Dish, transparent-red  
- 60x62x18 mm  
- [PS]

**Accessories**


910.66.1 Basket for dish washer for Inlay-Dishes



**Article No.** 66.Spez

**Description** Singe-Use-Inlay-Dish



**Details** - for MediDispenser No. 61/69  
- 60x62x18 mm  
- [PET]



**Article No.** 67

**Description** Inlay-Blister with 4 compartments

**Details** - for MediDispenser No. 61/69  
- 219x63x19 mm  
- [PET]




**Article No.** 341/349

**Description** MediDispenser 341

**Details** - MediDispenser No. 341  
with Sliding-Lid No. 49-69-349  
- 3 Divisions No. 342,  
divisible in max. 8 compartments  
- also available as imprinted Week-Set  
- 228x81x20 mm

**Accessories**

510.44xx Label with week days  
510.46 Label, short strip  
991 Week-Tray 41-61-341  
193TQ-41-61-341.1bl WeekBox 41-61-341, blue design  
191 Week-Container 41-61-341



## MediDispenser | the «large» ones


**Article No.** 141/149

**Description** MediDispenser 141

- Details**
- MediDispenser No. 141 with Sliding-Lid No. 149
  - 3 Divisions No. 142, divisible in max. 6 compartments
  - 228x78x23 mm

**Accessories**

- 510.14-8xx Label with dispensing periods
- 510.54 Label with dispensing times
- 510.13 Label for patients name and medication
- 510.13Wxx Label for medication with week days
- 510.13A4 Label "No. 13" for imprinting with MS-Word form
- 491 Week-Tray 141




**Article No.** 151/159

**Description** MediDispenser 151

- Details**
- MediDispenser No. 151 with Sliding-Lid No. 159
  - 3 Divisions No. 152, divisible in max. 7 compartments
  - 221x102x45 mm

**Accessories**

- 510.14-8xx Label with dispensing periods
- 510.54 Label with dispensing times
- 510.18 Label, extra-small
- 891 Storage Tray for 8 MediDispensers 151



**Article No.** 171/179

**Description** MediDispenser 171

- Details**
- MediDispenser No. 171 with Sliding-Lid No. 179-189
  - 3 Divisions No. 172, divisible in max. 8 compartments
  - 260x98x23 mm

**Accessories**

- 510.44xx Label, transparent (week days)
- 510.48 Label, long strip
- 193TQ-171.2gn WeekBox 171, green design




**Article No.** 181/189

**Description** MediDispenser 181

- Details**
- MediDispenser No. 181 with Sliding-Lid No. 179-189
  - 3 Divisions No. 182, divisible in max. 8 compartments
  - 260x98x39 mm

**Accessories**

- 510.44xx Label, transparent (week days)
- 510.48 Label, long strip





**Article No.** 161/169


**Description** MediDispenser 161

**Details**

- MediDispenser No. 161 with Sliding-Lid No. 169-269
- 4 Divisions No. 162, divisible in max. 12 compartments
- 363x77x23 mm

**Accessories**

- 510.54 Label with dispensing times
- 510.13 Label for patients name and medication
- 510.13Wxx Label for medication with week days
- 510.13A4 Label "No. 13" for imprinting with MS-Word form
- 691 Week-Tray 161-261



**Article No.** 261/269


**Description** MediDispenser 261

**Details**

- MediDispenser No. 261 with Sliding-Lid No. 169-269
- 4 Divisions No. 262, divisible in max. 12 compartments
- 363x77x45 mm

**Accessories**

- 510.54 Label with dispensing times
- 510.13 Label for patients name and medication
- 510.13Wxx Label for medication with week days
- 510.13A4 Label "No. 13" for imprinting with MS-Word form
- 691 Week-Tray 161-261



**Article No.** 271/279


**Description** MediDispenser 271

**Details**

- MediDispenser No. 271 with Sliding-Lid No. 279-289
- 4 Divisions No. 272, divisible in max. 12 compartments
- 338x67x23 mm

**Accessories**

- 510.54 Label with dispensing times
- 510.20 Label for week days and patients name
- 1091 Storage Tray for 8 MediDispensers 271-281



**Article No.** 281/289


**Description** MediDispenser 281

**Details**

- MediDispenser No. 281 with Sliding-Lid No. 279-289
- 4 Divisions No. 282, divisible in max. 12 compartments
- 338x67x45 mm

**Accessories**

- 510.54 Label with dispensing times
- 510.20 Label for week days and patients name
- 1091 Storage Tray for 8 MediDispensers 271-281




## MediDispenser | Inscription

**Article No.** 510.13

**Description** Label for patient's name and medication


- Details**
- for MediDispenser No. 11, 141, 161, 261
  - bundle: 350 labels
  - 75x48 mm
  - adhesive, removable



**Article No.** 510.13Wxx

**Description** Label for medication with week days


- Details**
- with imprint: «MONDAY..., SUNDAY»
  - sheet for 1 Week (bundle for 50 weeks)
  - for MediDispenser No. 11, 141, 161, 261
  - 75x48 mm
  - adhesive, removable
  - Replace **xx** in the article number with the desired language code:  
510.13Wde = Deutsch, 510.13Wfr = Français,  
510.13Wit = Italiano, 510.13Wen = English,  
510.13Wnl = Nederlands



**Article No.** 510.13A4

**Description** Label "No. 13" for imprinting with MS-Word form


- Details**
- A4-sheet: 10 labels (box with 2500 labels)
  - for MediDispenser No. 11, 141, 161, 261
  - 75x48 mm
  - MS-Word form for imprinting:  
[www.wiegand.ch](http://www.wiegand.ch)
  - adhesive, removable



**Article No.** 15/15SD

**Description** Label holder

- Details**
- with Sliding-Lid No. 15SD
  - suitable for Label No. 510.13, No. 510.13Wxx and No. 510.13A4
  - to be plugged on a Distribution-Tray
  - 77x55 mm



**Article No.** 64.738

**Description** Dispenser Marker

- Details**
- to mark MediDispenser, MediCups, Week-Sets etc.
  - with integrated eraser
  - dish washer proof ink






**Article No.** 510.46

**Description** Label, short strip


- Details**
- sheet with 11 labels (bundle of 20 sheets)
  - for MediDispenser No. 41, 61, 341
  - 191x11 mm
  - adhesive, removable



**Article Nr** 510.64-7

**Description** Label for medication, insert

- Details**
- for MediDispenser No. 61/69
  - 50 sheets with 7 labels each
  - paper: US-Legal form
  - 206x46 mm



**Article No.** 510.48

**Description** Label, long strip

- Details**
- sheet with 10 labels (bundle with 20 sheets)
  - for MediDispenser No. 171, 181
  - 235x12 mm
  - adhesive, removable




**Article No.** 510.17xx

**Description** Label for medication

- Details**
- A4 sheet with 8 labels
  - to be placed on Distribution Trays or in Storage Tray No. 1191
  - MS-Word form: [www.wiegand.ch](http://www.wiegand.ch)
  - 140x42 mm
  - Replace xx in the article number with the desired language code:  
510.17de = Deutsch, 510.17en = English, 510.17fr = Français

**Accessories**

17SD Label holder for label 510.17xx



**Article No.** 17SD

**Description** Label holder for label 510.17xx

- Details**
- to be placed on Distribution Trays or in Storage Tray No. 1191
  - 150x45 mm
  - [PET]

**Accessories**

510.17xx Label for medication



**Article No.** 510.44xx

**Description** Label with week days

- Details**
- with imprint: «MONDAY..., SUNDAY»
  - sheet of 2x7 labels, transparent
  - 41x8 mm
  - adhesive, removable
  - replace **xx** in the article number with the desired language code:  
510.44de = Deutsch, 510.44fr = Français,  
510.44it = Italiano, 510.44en = English,  
510.44nl = Netherlands



**Article No.** 510.14-8xx

**Description** Label with distribution periods

- Details**
- with imprint: «MORNING, NOON, EVENING, NIGHT, RESERVE»
  - sheet of 8x5 labels, transparent
  - 35x4 mm
  - adhesive, removable
  - replace **xx** in the article number with the desired language code:  
510.14-8de = Deutsch, 510.14-8fr = Français,  
510.14-8it = Italiano, 510.14-8en = English,  
510.14-8nl = Nederlands



**Article No.** 510.54

**Description** Label with distribution times


- Details**
- with imprint: «01:00, 02:00..., 24:00»
  - sheet of 35 labels, transparent
  - 13x13 mm
  - ideal for divisible MediDispenser
  - adhesive, removable



**Article No.** 510.18

**Description** Label extra-small


- Details**
- A4 sheets of 96 labels
  - for MediDispenser or MediCups
  - 40x10 mm
  - MS-Word form for imprinting:  
www.wiegand.ch
  - adhesive, removable



**Article No.** 510.20

**Description** Label small

- Details**
- A4 sheets of 39 labels
  - for MediDispenser No. 271, 281
  - 60x20 mm
  - MS-Word form for imprinting:  
www.wiegand.ch
  - adhesive, removable






**Article No.** 31

**Description** Distribution Tray simple

- Details**
- with 19 grooves
  - for MediDispensers  
No. 01/09, No. 11/19, No. 131/139,  
No. 41/49, No. 61/69 or No. 341/349
  - for 19 MediCups No. 71.x
  - 450x340 mm

**Accessories**  
16 Universal Mounting Plate  
for Distribution Trays

**Note** also available with  
enforced base, No. 31B




**Article No.** 51

**Description** Distribution Tray extended

- Details**
- with 17 grooves
  - for MediDispensers  
No. 01/09, No. 11/19, No. 131/139,  
No. 41/49, No. 61/69 or No. 341/349
  - for 35 MediCups No. 71.x
  - 450x340 mm

**Accessories**  
16 Universal Mounting Plate  
for Distribution Trays

**Note** also available with  
enforced base, No. 51B




**Article No.** 81

**Description** Distribution-Tray Universal,  
small

- Details**
- with 10 grooves
  - for all MediDispenser
  - for 50 MediCups No. 71.x
  - 300x400 mm

**Accessories**  
16 Universal Mounting Plate  
for Distribution Trays




**Article No.** 241

**Description** Distribution-Tray Universal,  
medium

- Details**
- with 16 grooves
  - for all MediDispenser
  - for 80 MediCups No. 71.x
  - 450x400mm

**Accessories**  
16 Universal Mounting Plate  
for Distribution Trays



**Article No.** 251


**Description** Distribution-Tray Universal, large

**Details**

- with 23 grooves
- for all MediDispenser
- for 115 MediCups No. 71.x
- 600x400 mm

**Accessories**

16 Universal Mounting Plate for Distribution Trays



### Trays | Storage Trays for MediDispensers

**Article No.** 891

**Description** Storage Tray for MediDispenser 151

**Details**

- for 8 MediDispensers No. 151/159
- 450x450x47 mm



**Article No.** 1091

**Description** Storage Tray for MediDispenser 271-281

**Details**

- for 8 MediDispensers No. 271/279 or No. 281/289
- 600x400x25 mm



**Article No.** 1191

**Description** Storage Tray for MediDispenser 131

**Details**

- for 8 MediDispensers No. 131/139 or 8x4 MediCups No. 71.x
- with plug-in for MediCups
- with 8 inscription grooves
- 400x204x45 mm

**Accessories**

510.17xx Label for medication  
 17SD Label holder for label 510.17xx  
 299-1199 Sliding-Lid for Trays 291 and 1191



### Week-Systems | Working aids

**Article No.** 295

**Description** Sliding-Lid-Opener

**Details**

- for Week-Tray No. 291 with MediDispenser No. 01/09
- consists of plug-in unit No. 296 and pusher No. 297

**Accessories**

01/09 MediDispenser 01  
 291 Week-Tray 01






## Week-Systems | Week-Trays with Sliding-Lids

**Article No.** 291

**Description** Week tray 01

**Details** - for 7 MediDispenser No. 01/09  
- 400x204 mm

**Accessories**  
299-1199 Sliding-Lid for Tray 291 and 1191  
510.46 Label, short strip




**Article No.** 591

**Description** Week-Tray 11

**Details** - for 7 MediDispenser No. 11/19  
- 400x245 mm

**Accessories**  
599 Sliding-Lid for Week-Tray 11  
510.46 Label, short strip



## Week-Systems | Week-Trays

**Article No.** 991

**Description** Week-Tray 41-61-341

**Details** - for 7 MediDispensers No. 41/49,  
No. 61/69 or No. 341/349  
- 600x260 mm

**Accessories**  
510.46 Label, short strip



**Article No.** 491

**Description** Week-Tray 141

**Details** - for 7 MediDispensers No. 141/149  
- 600x260 mm

**Accessories**  
510.46 Label, short strip



**Article No.** 691

**Description** Week-Tray 161-261

**Details** - for 7 MediDispensers No. 161/169  
or No. 261/269  
- 600x400 mm

**Accessories**  
510.46 Label, short strip



## Week-Systems | Week Container


**Article No.** 91

**Description** Week Container 11

**Details** - for 7 MediDispensers No. 11/19  
- 265x177x58 mm

**Accessories**

510.46 Label, short strip  
99 Sliding-Lid for Week-Container No. 91



**Article No.** 191

**Description** Week Container 41-61-341

**Details** - for 7 MediDispensers No. 41/49,  
No. 61/69 or No. 341/349  
- 265x177x87 mm

**Accessories**

510.46 Label, short strip  
199 Sliding-Lid for Week-Container No. 191



## Week-Systems | WeekBox

**Article No.** 193TQ-171.2gn

**Description** WeekBox 171, green design

**Details** - Week container for  
7 MediDispensers No. 171/179  
- foldable, washable  
- translucent white, green design  
- 267x160x102 mm

**Accessories**

510.95 Label for WeekBox




**Article No.** 193TQ-41-61-341.1bl

**Description** WeekBox 41-61-341, blue design

**Details** - Week container for 7 MediDispenser  
No. 41/49, No. 61/69 or No. 341/349  
- foldable, washable  
- translucent white, blue design  
- 138x84x229 mm

**Accessories**

510.95 Label for WeekBox  
7-41en/49 Week set 41 with week days  
7-61en/69 Week set 61 with week days  
7-341en/349 Week set 341 with week days



## MediDispenser Special | BlisterBox

**Article No.** 311

**Description** Dosys Dispenser

**Details** - Dispenser for blistered  
unit dose drugs  
for 7-14 days  
- 82x150x150 mm



**Article No.** 311T

**Description** BlisterBox Large

- Details**
- Dispenser for blistered unit dose drugs for 7-14 days
  - foldable, washable
  - translucent white
  - 78x162x139 mm



**Article No.** 312T

**Description** BlisterBox Compact

- Details**
- Dispenser for blistered unit dose drugs for 7 days
  - foldable, washable
  - translucent white
  - 78x128x107 mm



**Note** also available in white with imprint, No. 312W.1bl

**Article No.** 500.240XO.075

**Description** Holder for labels and blister bags

- Details**
- for blister bags
  - to be plugged on a Distribution-Tray
  - can be used as Label holder on SwissModul Trays and Baskets
  - 75x40 mm



## Storage | PharmaBox

**Article No.** 111/119

**Description** PharmaBox small

- Details**
- Container to store drugs in defined quantities
  - Sliding-Lid No. 119
  - 130x72x30 mm

**Accessories**

- 112e Label, mini-card
- 112R Mini-cards on rolls, white/yellow
- 112a Label holder for mini-cards
- 119weiss Sliding Lid, UV, for PharmaBox small
- H121-111 Support plate for PharmaBoxes small on shelves



**Article No.** 411/419

**Description** PharmaBox large

- Details**
- Container to store drugs in defined quantities
  - Sliding-Lid No. 419
  - 191x99x44 mm

**Accessories**

- 112e Label, mini-card
- 112R Mini-cards on rolls, white/yellow
- 112a Label holder for mini-cards
- 419w Sliding Lid, UV, for PharmaBox large
- H121-411 Support plate for PharmaBoxes large on shelves





## Tablet processing

**Article No.** 97

**Description** PillCrusher WIEGAND®

**Details**

- mortar for tablets
- 2-pieces, Ø 70 mm, height 105 mm
- with 1 MediCup, 25 ml

**Accessories**  
71.x MediCup



**Article No.** 96T

**Description** MediSplitter WIEGAND®

**Details**

- tablet cutter
- 85x35x25 mm
- blade: stainless steel
- with deposit compartment



**Article No.** 98T.1

**Description** MediCrusher WIEGAND®

**Details**

- mortar for tablets
- ergonomically designed
- with deposit compartment in the cover



**Article No.** 75.01

**Description** Dispenser Bag

**Article Info**

- for partitioned tablets or effervescent tabs
- 45x60 mm



## Drug storage | TRESTON

**Description** TRESTON Container, transparent

**Details**

- incl. Front Label and Label Protection
- ordering:
- No. 3010-1 = 92x300x82 mm
- No. 3015-1 = 132x300x100 mm
- No. 3020-1 = 186x300x82 mm
- No. 4010-1 = 92x400x82 mm
- No. 4015-1 = 132x400x100 mm
- No. 4020-1 = 186x400x82 mm

**Accessories**

- 10-D Division matching with No. 3010, No. 4010
- 15-D Division matching with No. 3015, No. 4015
- 20-D Division matching with No. 3020, No. 4020
- 10-E Front-Label matching with No. 3010, No. 4010
- 15-E Front-Label matching with No. 3015, No. 4015
- 20-E Front-Label matching with No. 3020, No. 4020
- 06-E Label for Divisions



## MediCups and accessories

**Article No.** 200.77.x

**Description** MediCups light 25 ml

- Details**
- medicine cup with graduation to 25ml, CE 0413
  - **x** to be ordered in different colours:  
 No. 200.77.1 = MediCup transparent-white  
 No. 200.77.2 = MediCup transparent-yellow  
 No. 200.77.3 = MediCup transparent-blue  
 No. 200.77.4 = MediCup transparent-red  
 No. 200.77.5 = MediCup transparent-green

**Accessories**  
200.78 Cover for MediCups



**Article No.** 71.x

**Description** MediCups 25 ml

- Details**
- medicine cup with graduation to 25ml, CE 0543
  - **x** to be ordered in different colours:  
 No. 71.1 = MediCup transparent-white  
 No. 71.2 = MediCup transparent-yellow  
 No. 71.3 = MediCup transparent-blue  
 No. 71.4 = MediCup transparent-red  
 No. 71.5 = MediCup transparent-green

**Accessories**  
72 Cover for MediCups  
910.71.1 Dish Washer Basket with insert  
64.746 Dispenser for MediCups



**Article No.** 64.746

**Description** Dispenser for MediCups

- Details**
- for 50 MediCups, 25 ml
  - to be fixed on the wall

**Accessories**  
71.x MediCups



**Article No.** 64.747

**Description** Dispenser for MediCup covers

- Details**
- Pic-Box Container, transparent, with Inox dispenser holder
  - to be fixed on the wall

**Accessories**  
72 Cover for MediCups



**Accessories for injections**

**Article No.** 64.679-2.0

**Description** Risk-Waste-Container large

- Details**
- E-SAFE, Universal-A-cover, yellow, 2.0 litre
  - with cannula-, scalpel and Insulin pen removal
  - flab against over-filling
  - closure for disposal
  - extra thick and secure material

**Accessories**

64.678 Holder for Risk-Waste-Container large for FlexLine Trolleys

Movix2.0 Movix, Vacuum-adhaesion-plate

E-TRAY45101 Syringe-Tray with adapter for small and large Risk-Waste-Containers

- Note**
- 0.5 litre container (No. 64.679-0.5) available
  - Vacuum-adhaesion-plate for 0.5 litre container (No. Movix0.5) available
  - Syringe-Tray with adapter for small Risk-Waste-Containers (E-TRAY45041)



**Article No.** E-Tray45101

**Description** Syringe-Tray

- Details**
- white, with Adapter for 0.5litre and for 2.0 litre Risk-Waste-Containers
  - with steril-unsteril-division
  - 350x250x40 mm




**Baskets for dish washer**

**Article No.** 910.71.1

**Description** Dish washer basket with MediCup insert


- Details**
- dish washer basket with lid
  - insert for 22 MediCups
  - 415x90x70 mm



**Article No.** 910.66.1

**Description** Dish washer basket with Inlay-Dish insert

- Details**
- dish washer basket with lid
  - insert for 18 Inlay-Dishes
  - 415x90x70 mm







### Transparent Sliding-Lids

Sliding-Lid production is a well guarded secret. The covers close the MediDispenser so firmly that no humidity can penetrate. A series of scientific tests in England showed outstanding results.

#### Material

PET: Polyethylene terephthalate is a polyester plastic (suitable for food). Contains no latex and no softeners. PET is UV-B, UV-C-FUV, UV-C-VUF, EUV and XUV resistant.

#### Recycling

- Emission free recycling
- Energy recycling (fuel) no reservations
- Very good ecological balance

#### Cleaning instructions

Washing	Clean with luke-warm water and mild detergent. (Hydrolysis)
Chemicals	Not allowed
Dishwasher	Not suitable
Sterilisation	Gas, radio-isotope (no steam!)

From approx. 45° C, WIEGAND® Sliding-Lids begin to open up their fold, due to the special manufacturing procedure. Thus the strength of the clamp reduces and a hermetic closure is no longer ensured.

### Blue MediDispenser and Trays

All blue WIEGAND® products are manufactured deliberately in elaborate deep-drawing procedures. Material formed in this way becomes very resistant, and even withstands dropping. The "blue line" is light sensitive. Avoid UV-exposure.

#### Material

Polystyrene (food foil). Contains no latex, fluorine, cadmium, lead, bromine or chlorine.

#### Recycling

- Emission free recycling
- Energy recycling (fuel) no reservations
- Very good ecological balance

#### Cleaning reference

Chemicals	with frequent use of caustic solutions, alcohol and acids, the plastic becomes brittle
Dishwasher	at 55° C (eco-program) without any problem, at 65° C maximum 45 min. running time at 95° C (industrial rinsing machines) 2 min. running time
Sterilisation	gas, radioisotopes, (no steam!)

WIEGAND® Dispensers have been certified by various European states.



**SUMINISTROS PUNTO  
GERIÁTRICO**



**TELEFONO 654 512487**